

SNOREMENDER® 3 INSTRUCTIONS

These instructions, and our FAQ, are available at www.snoremenders.com

When SnoreMender cannot be used

The SnoreMender can not be used in the following situations:

- If you have an excessively receding jaw (overbite 5mm (3/16 inch) or more) or a protruding lower jaw.
- If you suffer from Gum Disease / Jaw Pain / TMJ dysfunction (TMD) / Loose Teeth.
- If you have Dentures / Implants / Braces / Bridgework (not including crowns) / Other Weaknesses in your front 6 upper or front 8 lower teeth you should check with us or your dentist, to make sure they are strong enough to support the SnoreMender.
- People under 18 years old MUST NOT use the SnoreMender.

How to Use

If you experience prolonged discomfort, bleeding, loose or moving teeth, jaw pain, bite changes or other problems related to using your SnoreMender, please immediately discontinue use and contact us for advice.

With proper use there is no permanent effect of the SnoreMender on jaw or teeth - please make sure you read the usage instructions carefully.

To Use: Before going to sleep, fold the SnoreMender and place in the mouth.



← Fits over Upper teeth

← Fits behind Lower teeth

The SnoreMender may hold your mouth open wide at first. This goes away in a couple of days, but you can speed it up by putting some boiling water in a cup and holding the hinges in the water for 30 seconds with it folded closed.

The SnoreMender works by moving the lower jaw slightly forwards which can be felt as a forward pressure on the lower teeth, indicating correct use.

The SnoreMender should work immediately, but the jaw and teeth may ache for the first few weeks, which is normal as your body adjusts. This usually goes away within

a few minutes of removing the SnoreMender but teeth may feel sensitive for longer at first.

The key is to persevere - give yourself time to get used to it. Our staff includes long time users of SnoreMenders – if you have any difficulties, please see our FAQ at www.snoremenders.com or email us at help@snoremenders.com

If the teeth continue to ache significantly, after the first 8 weeks of use, please check the "If the SnoreMender is too strong" section of the instructions.

Occasionally at night your SnoreMender may fall out. This is normal and should happen less over time - if it keeps falling out, please see adjustments below.

Sometimes the SnoreMender may cause excess saliva or a dry mouth - your mouth will usually adjust within two to four weeks.

Cleaning and Maintaining SnoreMender

Warning - never overheat your SnoreMender – Max 100°C / 212°F !!!

1. Wash by hand in warm water, with soap, and rinse.
2. If you want to keep more springiness in the hinges, once clean, hold the SnoreMender opened out flat and run hot tap water over the hinges for 30 seconds to soften them slightly. Then dry and store in the box.

If the SnoreMender is rubbing the gums

You can use a nail file and file the edges down, or nail scissors to trim them. You will need to press firmly with the file. It is best to do one half first, so you can see how much you have filed. About 1 mm (1/32 inch) is usually enough.

You can smooth any rough edges by using the nail file gently.

If it rubs elsewhere, please email help@snoremenders.com for advice.

If the SnoreMender Slips out of place

Sometimes the lower part of the SnoreMender may be a little too wide, so it will not fit behind the lower teeth without squeezing it into place, and it will then sometimes 'pop off the teeth'. In this case, the following adjustment should help:

1. Put 3 cups of **boiling** water in a bowl.
2. Hold the lower front half of the device in the water for 90 seconds.
3. Squeeze the lower front part together firmly to make it narrower.
4. Hold in position and put in cold water for 30 seconds to cool.


If you make it too narrow, you can just put it back in the hot water for a few seconds and it will return to the original shape.

If the slippage is elsewhere, please email help@snoremenders.com for advice.

If the SnoreMender Is Too Strong

There is always an adjustment period when you start using an anti-snoring mouthpiece so, if possible, we recommend you use it without adjustment for 1-2 weeks. However, if it is too strong to wear, or it is still making the teeth ache after two weeks use, the following adjustment will reduce the strength.

The goal is to lengthen the middle of the upper tray (typically by 5mm / a sixth of an inch), to reduce the jaw advancement and pressure on the teeth and jaw.




1. Put 3 cups of **boiling** water in a bowl.
2. Put the upper front part in the water for 60 seconds.
3. Put your fingers inside the top half of the SnoreMender and pull it apart, to stretch the front area of the upper part.
4. Hold in position for 60 seconds while it cools.
5. If you stretch it too much you can put the stretched part back in the hot water and it should shrink back to the original size.

If that does not help, please email help@snoremenders.com for advice.

If the SnoreMender is Not Working

If the SnoreMender is not stopping snoring it may not be extending the lower jaw enough to open the airway. This adjustment will increase jaw advancement:

1. Put 3 cups of **boiling** water in a bowl.
2. Hold the lower front piece (NOT the hinges) in the water for 60 seconds, while it softens.
3. Pull the lower tray outwards (by putting thumb through middle and pulling away from hinges) and hold for 60 seconds. You should feel it lengthening.
4. You should see it is now narrower and elongated.
5. You can try it several times, until you get the length required - start cautiously and get the feel of how hard to pull it.



If that does not help, please email help@snoremenders.com for advice. We make an extra strength version of the SnoreMender and there is no charge to switch to it if required.

Other fitting problems - general moulding

Sometimes the teeth may be misaligned, or there may be parts of the mouth that the SnoreMender is rubbing. The following adjustment will make the SnoreMender fit better:

1. Put 3 cups of **boiling** water in a bowl.
2. Put the whole device in the water for 60 seconds.
3. Put the device in the mouth and bite down gently, to force it into place.
4. Hold in position for 90 seconds while it cools.

If that does not help, please email help@snoremenders.com for advice.

Overactive Gag Reflex - Acclimatisation

The following steps should help while you acclimatise:

- Warm it slightly in warm water before wearing.
- Hold it in place with your fingers, so your body knows you are in control of it.
- Start a few seconds at a time and build up your ability to not gag.

Replacing Your SnoreMender

Eventually your SnoreMender will need replacing - typical lifetimes are between 9 and 18 months. If you notice it becoming less effective, please visit www.snoremenders.com to purchase a replacement.